

培养基配方表 (Media Formulation)

名称: MEM 必需氨基酸溶液 50X

货号: S230JV

组分 (Components)	分子量 (MW)	浓度 (mg/L)	摩尔浓度 (mM)
无机盐 (Inorganic Salts)			
L-Arginine	174.2	6320	36.28013777
L-Cystine	240.3	1200	4.993757803
L-Histidine monohydrochloride hydrate	209.63	2100	10.01765015
L-Isoleucine	131.17	2620	19.97407944
L-Leucine	131.17	2620	19.97407944
L-Lysine hydrochloride	182.65	3625	19.84670134
L-Methionine	149.21	755	5.059982575
L-Phenylalanine	165.19	1650	9.988498093
L-Threonine	119.12	2380	19.97985225
L-Tryptophan	204.23	510	2.497184547
L-Tyrosine	181.19	1800	9.934323086
L-Valine	117.15	2340	19.97439181